





Cover the lid, rotate the lid clockwise in parallel, until the clasps slide into the slots to lock.




Rotate the lid counterclockwise in parallel until the clasps slip out of the slots to unlock, then the lid can be taken away.

## Product Range

| Cat. No. | Product Description           | /Pack | /Case |
|----------|-------------------------------|-------|-------|
| 752101   | 90×15 mm Petri Dish with Lock | 20    | 500   |
| 722311   | 65×15 mm Petri Dish with Lock | 20    | 500   |

## Technical Drawing

## 90×15 mm Petri Dish with Lock:




## Head Office

Email: [info@cell-nest.com](mailto:info@cell-nest.com)  
 Online: [www.cell-nest.com](http://www.cell-nest.com)

## Oversea

NEST USA (New Jersey/ Phoenix)  
 NEST scientific 株式会社 (Yokohama, Japan)  
 NEST Scientific Europe B.V (Netherlands)  
 Nest Scientific (MENA) FZE (Sharjah , United Arab Emirates)

## 65×15 mm Petri Dish with Lock:



### Head Office

Email: [info@cell-nest.com](mailto:info@cell-nest.com)  
Online: [www.cell-nest.com](http://www.cell-nest.com)

### Oversea

NEST USA (New Jersey/ Phoenix)  
NEST scientific 株式会社 (Yokohama, Japan)  
NEST Scientific Europe B.V (Netherlands)  
Nest Scientific (MENA) FZE (Sharjah , United Arab Emirates)

## Product picture



### Head Office

No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi,  
Jiangsu, China  
Tel : +86+ 510-6800 6788 Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
Online : [www.cell-nest.com](http://www.cell-nest.com)

### Oversea

NEST USA (New Jersey/ Phoenix)  
NEST scientific 株式会社 (Yokohama, Japan)  
NEST Scientific Europe B.V (Netherlands)  
Nest Scientific (MENA) FZE (Sharjah , United Arab  
Emirates)