

Product Name:	Biofactory
Product Application:	Nest BioFactory™ Systems are compact, multi-layer, single-use cell culture systems designed for easy scale-up cell culture applications, can be used for medium to large scale research and industrial-scale production such as vaccines, monoclonal antibodies and biopharmaceuticals. It is suitable for adherent cell culture. Linear amplification will not change the kinetic conditions of cell growth.
Product Materials:	Polystyrene (PS), meets USP, Class VI.
product Specifications:	Temperature range: stored at room temperature Shelf life: 3 Years after date of production (ensure package is in good)
Sterilization:	Yes. Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
Pyrogens:	Non-Pyrogenic
RNase/DNase:	DNase/RNase free
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Particle Test: qualified Hydrophilic test: qualified Cell culture test: qualified Packaging inspection: qualified
Features	1. Made of high clarity, Strictly selected polystyrene conforming to USP Class VI standard. 2. Class 10,000 clean production workshop, process has been repeatedly verified, strict quality management, to ensure stable product quality. 3. The irradiation sterilization process has passed ISO 13485 and ISO 11137 quality system certification. 4. TC treated, good for adherent cell attachment and growth. 5. Ultrasonic welding without extrinsic ingredients, without any chemical addition, eliminating the risk of unidentified leachables. Excellent welding line design can avoid producing particles during welding and increases the mechanical strength. 6. Vent caps with 0.22um hydrophobic filters for gas exchange without contamination. 7. Wide mouth is applicable to pour culture medium directly. Narrow mouth is applicable to operate with the aseptic pipeline. 8. Growth kinetics are the same as cells grow in the cell culture flasks. Applicable to large-scale culture adherent cells. 9. Large growing area. Just a single operate, you can cultivate a large number of cells and

Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

reduce contamination risk.

10. Biosafety testing Refer to <<State Food and Drug Administration National Standards for Packaging Materials and Containers in Direct Contact with Drugs (Series 6)>>

11. Applicable to automatic machine.

12. Standardize operations to reduce batch differences.

13. Double-layer sterile bag, clear lot number for batch traceability.

Product Range

Cat.No.	Layer	Cap Style / Case	Cultivation area (cm ²)	/ Case
771001	1	16 sterile wide vent caps (with 16 sterile wide plug seal caps included separately)	632	8
771101	2	16 sterile wide vent caps (with 16 sterile wide plug seal caps included separately)	1264	8
771704	4	8 sterile wide vent caps (with 8 sterile wide plug seal caps included separately)	2528	4
771204	5	8 sterile wide vent caps (with 8 sterile wide plug seal caps included separately)	3160	4
771604	6	8 sterile wide vent caps (with 8 sterile wide plug seal caps included separately)	3792	4
771302	10	12 sterile wide vent caps (with 12 sterile wide plug seal caps included separately)	6320	6
771503	15	4 sterile wide vent caps (with 4 sterile wide plug seal caps included separately)	9480	2
771403	40	4 sterile wide vent caps (with 4 sterile wide plug seal caps included separately)	25280	2
772001	1	16 sterile narrow vent caps (with 16 sterile narrow plug seal caps included separately)	632	8
772101	2	16 sterile narrow vent caps (with 16 sterile narrow plug seal caps included separately)	1264	8
772704	4	8 sterile narrow vent caps (with 8 sterile narrow plug seal caps included separately)	2528	4
772204	5	8 sterile narrow vent caps (with 8 sterile narrow plug seal caps included separately)	3160	4
772302	10	12 sterile narrow vent caps	6320	6

Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

		(with 12 sterile narrow plug seal caps included separately)		
772403	40	4 sterile narrow vent caps (with 4 sterile narrow plug seal caps included separately)	25280	2
773001	1	8 sterile wide vent caps & 8 sterile narrow vent caps (with 8 sterile wide plug seal caps & 8 narrow plug seal cap included separately)	632	8
773101	2	8 sterile wide vent caps & 8 sterile narrow vent caps (with 8 sterile wide plug seal caps & 8 narrow plug seal cap included separately)	1264	8
773204	5	4 sterile wide vent caps & 4 sterile narrow vent caps (with 4 sterile wide plug seal caps & 4 narrow plug seal cap included separately)	3160	4
773302	10	6 sterile wide vent caps & 6 sterile narrow vent caps (with 6 sterile wide plug seal caps & 6 narrow plug seal cap included separately)	6320	6
773403	40	2 sterile wide vent caps & 2 sterile narrow vent caps (with 2 sterile wide plug seal caps & 2 narrow plug seal cap included separately)	25280	2

Head Office


Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea


NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

Technical Drawing


Wide mouth:



Narrow mouth:



2 Wide Caps:



Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)


1 Layer




2 Layer




4 Layer




5 Layer




6 Layer



10 Layer



15 Layer




Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com


Oversea

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

40 Layer



2 Narrow Cap:



Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com


Oversea

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)


1 Layer




2 Layer




4 Layer




5 Layer



10 Layer



40 Layer




Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

1 Wide Cap & 1 Narrow Cap :




1 Layer




2 Layer



5 Layer



10 Layer




Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea


NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

40 Layer



Degree Scale

1 Layer




Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com


Oversea

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)


2 Layer



5 Layer



10 Layer




Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com


Oversea

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)


40 Layer



Wide Vent/ Plug Seal Caps



Narrow Vent/ Plug Seal Caps



Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)