




Product Name:	250 mL & 500mL Centrifuge Tubes	
Product Application:	The 250mL/500mL centrifuge tubes can process a large number of samples at one time and are suitable for separation and harvest of large numbers of cells, supernatants, bacteria, yeast and tissue samples, therefore they are widely used in the field of life sciences and clinical experiments and can meet the requirements of biological analysis. 1. Separation of sediments after chemical reactions for harvest and application. 2. Harvest of biomacromolecules, inorganics and organics. 3. In biochemistry and other biological fields, usually used to centrifuge and collect a large volume of cells, biomacromolecules and supernatants.	
Product Materials:	PP (tube) + HDPE (cap), Meets USP, Class VI standards	
product Specifications:	Temperature range: stored at room temperature Shelf life: 5 Years after date of production (ensure package is in good)	
Sterilization:	Yes. Sterilized by E-beam, Sterility Assurance Level: SAL=10⁻⁶. The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137	
Pyrogens:	Non-Pyrogenic	
RNase/DNase:	DNase/RNase free	
BSE/TSE:	These products are deemed animal free.	
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Centrifugation test: qualified Freeze test: qualified Negative pressure test: qualified Packaging inspection: qualified	
Features	1. Unique structure design of the cap provides double security, further improve the seal ability. 2. Increased the thickness of conical bottom in order to improve the maximum centrifugal limit.   3. Max. centrifugal limit 7000xg. 4. Tolerable temperature range: -80°C~120°C.	

5. Clear scale for easy viewing.
6. DNase, RNase and pyrogen free.
7. Electron beam sterilization, SAL=10⁻⁶.
8. Tubes and cap free of heavy metal.


Product Range

Cat. No.	Product Description	/Pack	/Case
622001	250mL Centrifuge Tube, Bulk, Sterilized	6	60
623001	500mL Centrifuge Tube, Bulk, Sterilized	6	36
622002	250mL Centrifuge Tube, PS Rack, Sterilized	6	24
623002	500mL Centrifuge Tube, PS Rack, Sterilized	4	24

Technical Drawing



500ML



250ML


Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

500ML



250ML



Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)